



LBMA ASSAYING & REFINING Conference 2023



**Improving the sensitivity, precision and accuracy of
ICP-OES analyses to fulfill standardization regulations**

Dr. Dirk Wüstkamp, Senior Product Specialist ICP-OES, SPECTRO Analytical Instruments GmbH

ISO/TC 174 - Jewellery and precious metals

Created 1978

Secretary: DIN Germany

Chairperson: Ph.D. Jonathan J. Jodry, Switzerland

ISO/TC 174/WG 1 Methods for determining fineness

ISO/TC 174/WG 2 Diamonds, gemstones and related products

18 Participating members (countries)

23 Observing members (countries)

20 published ISO standards, under the direct responsibility of ISO/TC 174

8 ISO standards under development, under the direct responsibility of ISO/TC 174

ISO/TC 174/WG 1 Methods for determining fineness

ISO 9202 –Fineness of precious metal alloys

Precious metal	Fineness ‰	Recommended method
Ag	50 - 990	ISO NP 19919 (DIN 32562)
	333 - 990	ISO 11427 + 13756
	800 - 990	ISO 11427 + 13756
	≥ 999	ISO 15096
	999.9	ISO 15096
	999.99	ISO 15096
Au	100 - 999	ISO 11426
	≥ 999	ISO 15093
	999.9	ISO 15093
	999.99	ISO 15093
	999.999	under development (ICP-MS and GD-MS)
Pt	50 - 990	ISO 11210
	≤ 990	ISO 11494
	≥ 999	ISO 15093
Pd	50 - 990	ISO 11490
	≤ 990	ISO 11495
	≥ 999	ISO 15093


Dual Side-On Interface – Principle

- ▀ Since the light emitted into two directions is transferred into the optical system the sensitivity is in average improved by a factor of two
- ▀ In addition, the vertical torch provides high stability and freedom from matrix effects since the affected zones in the plasma are blanked out

- ▀ Pros and Cons:

- **High stability**
- **High matrix tolerance**
- **High linear dynamic range**
- **Freedom from matrix effects**
- **All with ONE measurement**

- Only a somewhat lower (factor 2) sensitivity compared to a direct light path axial plasma observation



ISO 15093 - DSOI - Limit of detection in 2 % Gold

Element	LOD in solution mg/l	LOD solid material mg/kg	Element	LOD in solution mg/l	LOD solid material mg/kg
Ag 328.068	0,014	0,70 (10)	Mo 202.095	0,0025	0,13
Al 167.078	0,0002	0,01	Ni 221.648	0,0007	0,04 (3)
As 193.759	0,013	0,65 (30)	P 177.495	0,0017	0,09
B 182.641	0,0022	0,11	Pb 168.215	0,0018	0,09 (10)
Ba 455.404	0,0003	0,02	Pd 324.270	0,013	0,65 (10)
Be 313.042	0,00005	0,003	Pt 177.708	0,012	0,6
Bi 223.061	0,013	0,65 (10)	Rh 343.489	0,033	1,65
Ca 396.847	0,001	0,05	Ru 240.272	0,005	0,25
Cd 228.802	0,0005	0,03	S 182.034	0,0095	0,48
Co 230.786	0,0001	0,005	Sb 206.833	0,0054	0,27
Cr 283.563	0,0013	0,07 (3)	Se 196.090	0,0047	0,235
Cu 324.754	0,0026	0,13 (10)	Sn 189.991	0,0036	0,18 (10)
			Si 251.612	0,0054	0,27 (10)
Fe 239.562	0,001	0,05 (10)	Sr 421.552	0,0002	0,009
Ge 164.919	0,003	0,15	Ti 334.941	0,00037	0,02
Hg 194.227	0,0039	0,20	Tl 190.864	0,009	0,45
Ir 183.250	0,0038	0,19	V 292.402	0,001	0,05
Li 670.780	0,004	0,2	W 239.709	0,03	1,5
Mg 280.270	0,00034	0,02 (10)	Zn 213.856	0,0003	0,02
Mn 259.373	0,00022	0,01 (3)	Zr 343.823	0,0037	0,19

What purity level can be assessed by ICP-OES ?

$$\sum LOD < 10 \text{ mg/kg}$$

99.999



ASTM B 562 - 95

$$\sum LOD < 4 \text{ mg/kg}$$

99.9996

*(5) maximum Limits ASTM 99.995

ISO 15096 – DSOI - Limit of detection in 5 % Silver

Element	LOD in solution	LOD solid material	Element	LOD in solution	LOD solid material
	mg/l	mg/kg		mg/l	mg/kg
Al 167,078	0,0006	0,01	Mo 202,095	0,0007	0,01
As 189.042	0,003	0,07	Nb 316,340	0,002	0,05
Au 174.050	0,008	0,15	Ni 231.604	0,005	0,09
B 249.773	0,0013	0,03	P 177.495	0,002	0,04
Ba 233.527	0,0003	0,006	Pb 168.215	0,002	0,04 (10)
Be 234.861	0,000006	0,0001	Pd 324.270	0,025	0,50 (10)
Bi 223.061	0,005	0,09 (5)	Pt 214.423	0,002	0,05
Ca 396.847	0,002	0,05	Rh 233.477	0,015	0,30
Cd 214.438	0,0003	0,005	*Ru 240.272	0,024	0,47
Co 238.892	0,006	0,12	S 182.034	0,005	0,09
Cr 267.716	0,001	0,02	Sb 217.581	0,004	0,08
Cu 219.226	0,004	0,08 (100)	Se 196.090	0,004	0,07 (5)
Fe 259.941	0,0006	0,01 (10)	Sn 189.991	0,001	0,03
Ga 141.444	0,007	0,14	Sr 407.771	0,00008	0,002
Ge 164.919	0,004	0,08	Ta 226.230	0,004	0,09
Hf 264,141	0,002	0,04	Te 214.281	0,004	0,08 (5)
Hg 184.950	0,002	0,03	Ti 334,187	0,001	0,03
In 158.637	0,005	0,09	Tl 190.864	0,004	0,09
*Ir 183.250	0,002	0,05	V 292.402	0,001	0,03
Li 670.780	0,001	0,03	W 220.448	0,005	0,10
Mg 279.553	0,0004	0,008	Zn 213.856	0,0006	0,01
Mn 260.569	0,0003	0,005	Zr 339.198	0,002	0,04

What purity level can be assessed by ICP-OES ?

$$\sum LOD < 4 \text{ mg/kg}$$

99.9996



ASTM B413-97a


$$\sum LOD < 1 \text{ mg/kg}$$

99.9999

* (5) maximum Limits ASTM 99.99

ISO 15096 - DSOI – Spike Recovery in 5 % Silver

Spike Recovery [%]




Spike Recovery Multi Element Standards

- **CCS 2** 1,0mg/l \triangleq 20 mg/kg
- **CCS 5** 0,25mg/l \triangleq 5 mg/kg
- **Merck VI** 0,025mg/l \triangleq 0,5 mg/kg



Results – round robin ISO 15096 - Ag

\bar{y}			
Average	s	z-score with	z-score
[%]	[%]	outliers	without
			outliers
999.662	0.005	-1.62	-1.64
999.664	0.004	-1.56	-1.57
999.688	0.003	-0.60	-0.59
999.691	0.003	-0.48	-0.46
999.699	0.009	-0.15	-0.12
999.703	0.001	-0.01	0.02
999.703	0.005	0.03	0.06
999.710	0.005	0.31	0.35
999.711	0.031	0.34	0.37
999.717	0.007	0.56	0.61
999.717	0.003	0.57	0.61
999.735	0.013	1.30	1.36
999.736	0.003	1.32	1.38
Great average		999.703	999.702



ISO 15093 - DSOI - Limit of detection in 2 % Palladium

Element	LOD in solution	LOD solid material	Element	LOD in solution	LOD solid material
	mg/l	mg/kg		mg/l	mg/kg
Ag 328.068	0,003	0,13 (100)	Mo 281.615	0,003	0,16
Al 167.078	0,003	0,16 (50)	Nb 269.706	0,003	0,14
As 228.812	0,014	0,68	Ni 221.648	0,001	0,06 (50)
Au 201.265	0,006	0,28 (100)	P 213.618	0,010	0,48
B 249.773	0,003	0,15	Pb 168.215	0,042	2,10 (50)
Ba 455.404	0,0004	0,02	Pt 177.708	0,004	0,21
Be 313.042	0,00001	0,0003	Re 227.525	0,004	0,19
Bi 223.061	0,013	0,64	Rh 343.489	0,025	1,27
Ca 393.366	0,001	0,04 (50)	Ru 240.272	0,007	0,37
Cd 214.438	0,001	0,03	S 182.034	0,006	0,29
Co 230.786	0,002	0,08 (10)	Sb 217.581	0,010	0,49 (20)
Cr 267.716	0,002	0,10 (10)	Se 196.090	0,063	3,13
Cu 327.396	0,004	0,21 (50)	Si 251.612	0,012	0,58 (50)
Fe 261.187	0,004	0,18 (50)	Sn 147.516	0,058	2,92 (50)
Ga 141.444	0,021	1,03	Sr 407.771	0,0001	0,01
Ge 209.426	0,011	0,55	Ta 240.063	0,011	0,55
Hf 277.336	0,004	0,18	Te 214.281	0,010	0,49
Hg 194.227	0,002	0,11	Ti 336.121	0,001	0,06
In 158.637	0,017	0,87	Tl 132.171	0,070	3,49
Ir 183.250	0,005	0,25	V 292.402	0,001	0,06
Li 670.780	0,003	0,13	W 207.911	0,004	0,21
Mg 279.553	0,002	0,08 (50)	Zn 213.856	0,0004	0,02 (25)
Mn 260.569	0,0004	0,02 (10)	Zr 343.823	0,002	0,08

What purity level can be assessed by ICP-OES ?

$$\sum LOD < 23 \text{ mg/kg}$$

99.997



ASTM B589-94


$$\sum LOD < 10 \text{ mg/kg}$$

99.999

*(25) maximum Limits ASTM 999.5

ISO 15093 - DSOI – Spike Recovery in 2 % Platinum

Spike Recovery [%]



Spike Recovery Multi Element Standards

- CCS 2 0,025mg/l \triangleq 1,25 mg/kg
- CCS 5 0,025mg/l \triangleq 1,25 mg/kg
- Merck VI 0,025mg/l \triangleq 1,25 mg/kg



ISO 15093 - DSOI - Limit of detection in 2 % Platinum

Element	LOD in solution mg/l	LOD solid material mg/kg	Element	LOD in solution mg/l	LOD solid material mg/kg
Ag 338.289	0,014	0,72 (30)	Mo 281.615	0,004	0,19 (40)
Al 396.152	0,016	0,78 (40)	Nb 316.340	0,004	0,18
As 197.262	0,034	1,72 (20)	Ni 231.604	0,002	0,11 (10)
Au 267.595	0,007	0,34 (50)	P 213.618	0,018	0,90
B 249.677	0,005	0,26	Pb 220.353	0,015	0,74 (10)
Ba 455.404	0,0004	0,02	Pd 340.458	0,017	0,87 (50)
Be 313.042	0,000004	0,0002	Rh 343.489	0,040	2,01 (50)
Bi 223.061	0,017	0,87 (20)	Ru 245.644	0,012	0,59 (20)
Ca 396.847	0,001	0,04 (30)	S 182.034	0,011	0,57
Cd 226.502	0,001	0,04	Sb 217.581	0,030	1,50 (20)
Co 238.892	0,002	0,10	Se 196.090	0,012	0,60
Cr 283.563	0,001	0,06 (10)	Sn 147.516	0,012	0,61 (20)
Cu 327.396	0,004	0,21 (40)	Sr 421.552	0,0002	0,01
Fe 261.187	0,004	0,22 (50)	Ta 240.063	0,006	0,31
Ga 417.206	0,032	1,61	Te 214.281	0,030	1,52 (40)
Ge 164.919	0,011	0,56	Ti 334.941	0,001	0,03
Hg 194.227	0,009	0,46	Tl 132.171	0,070	3,48
Ir 183.250	0,009	0,47 (50)	V 292.402	0,001	0,07
Li 670.780	0,002	0,09	W 220.448	0,009	0,46
Mg 280.270	0,001	0,07 (30)	Zn 213.856	0,001	0,07 (20)
Mn 257.611	0,0004	0,02 (10)	Zr 339.198	0,001	0,07

What purity level can be assessed by ICP-OES ?

$$\sum LOD < 24 \text{ mg/kg}$$

99.997



ASTM B561-94


$$\sum LOD < 14 \text{ mg/kg}$$

99.998

*(25) maximum Limits ASTM 999.9

ISO 15093 - DSOI – Spike Recovery in 2 % Platinum

Spike Recovery [%]



Spike Recovery Multi Element Standards

- CCS 2 0,025mg/l \triangleq 1,25 mg/kg
- CCS 5 0,025mg/l \triangleq 1,25 mg/kg
- Merck VI 0,025mg/l \triangleq 1,25 mg/kg



ARCOS EOP - Limit of detection in 5 % Gold

Element	LOD in solution	LOD solid material	Element	LOD in solution	LOD solid material
	mg/l	mg/kg		mg/l	mg/kg
Ag 328.068	0,0013	0,026	Mn 259.373	0,0003	0,006
Al 167.078	0,0006	0,01	Mo 202.095	0,0033	0,07
As 193.759	0,037	0,74	Ni 221.648	0,0021	0,04
B 182.641	0,0023	0,05	Pb 168.215	0,0104	0,21
Ba 455.404	0,0003	0,006	Pd 324.270	0,014	0,28
Bi 223.061	0,0113	0,23	Pt 177.708	0,0086	0,17
Cd 228.802	0,0008	0,02	Re 221.426	0,0041	0,08
Ce 413.765	0,017	0,34	Rh 343.489	0,0234	0,47
Co 230.786	0,0022	0,04	Ru 240.272	0,0072	0,14
Cr 283.563	0,0011	0,02	Sb 206.833	0,0102	0,20
Cu 324.754	0,0011	0,02	Sc 335.373	0,0007	0,01
Fe 239.562	0,0014	0,03	Se 196.090	0,0182	0,36
Ga 141.444	0,0111	0,22	Sn 189.991	0,0094	0,19
Ge 164.919	0,0082	0,16	Sr 407.771	0,0001	0,002
Hf 232.247	0,0036	0,07	Ta 268.517	0,012	0,24
Hg 184.950	0,0147	0,29	Te 214.281	0,031	0,62
In 158.637	0,0082	0,16	Ti 334.941	0,0006	0,01
Ir 183.250	0,0225	0,45	Tl 190.864	0,0204	0,41
Li 670.780	0,0034	0,07	Zn 206.200	0,0015	0,03
Mg 280.270	0,0001	0,002	Zr 343.823	0,0021	0,04

ISO NWIP : 99.9999 possible by ICP-OES ?

What purity level can be assessed by ICP-OES ?

$$\sum LOD < 7 \text{ mg/kg}$$

99.9993



ASTM B 562 - 95

$$\sum LOD < 2 \text{ mg/kg}$$

99.9998


Bracketing - Highest precision and accuracy

Procedure described in :

- **ISO 11494** **Platinum in alloys**
- **ISO 11495** **Palladium in alloys**
- **ISO NP 19919 (DIN 32562)** **Silver in alloys**
- **NWIP** **Gold in alloys**

General requirements:

- Mandatory use of an internal standard.
- High precision, weight-based preparations (analyte and internal standard).
- Repeated measurements (“cycles”) of the sample and two standards of a slightly lower and higher concentration.





BRACKETING



BRACKETING

ISO 11494 - ARCOS DSOI - Bracketing


High precision for main components

Sample	Type	Pt 306.471	
		mg/l	
	1	173,535	
	2	173,559	
	3	173,597	
	4	173,565	
	5	173,541	
RRT Platinum (liquid)	<x>	173,559	
	sd	0,024	
	rsd	0,014	
<i>EW (g)</i>	<i>0,09128</i>		
<i>Vol (ml)</i>	<i>500</i>		
		Pt 306.471	
		%	
RRT Platinum (solid)		95,070	
	RSD(%)	0,014	



Consideration of the expanded uncertainty of measurement ISO 11494

Consideration of the expanded uncertainty of measurement - ISO 11494



Expanded Uncertainty U_{cs}
 $U_{CS} = 2 * U_c = 0,10 \%$
 U_c = combined uncertainty

Multi-element bracketing for alloys

Sample	Typ	Pt %	Pd %	Cu %	Total %
Sample 1-1	1	60.33	10.10	29.57	99.99
	2	60.32	10.09	29.58	100.00
	3	60.29	10.08	29.59	99.96
	4	60.31	10.10	29.59	99.99
	5	60.28	10.08	29.61	99.97
	<x>	60.31	10.09	29.59	99.98
	rsd	0.033	0.102	0.072	0.046
Sample 1-2	1	60.312	10.083	29.61	100.01
	2	60.31	10.108	29.58	100.00
	3	60.311	10.082	29.59	99.98
	4	60.32	10.101	29.59	100.01
	5	60.334	10.082	29.59	100.00
	<x>	60.32	10.09	29.59	100.00
	rsd	0.017	0.124	0.018	0.012
Sample 1-3	1	60.316	10.09	29.56	99.97
	2	60.331	10.125	29.55	100.01
	3	60.292	10.091	29.59	99.97
	4	60.294	10.111	29.57	99.98
	5	60.266	10.09	29.63	99.99
	<x>	60.30	10.10	29.58	99.98
	rsd	0.05	0.155	0.088	0.022



RRT for 3 component alloy

Same procedure as in single element bracketing

	Pt %	Pd %	Cu %	Total %
Lab 1	60.31	10.09	29.59	99.99
Lab 2	60.38	10.04	29.58	100.00
Lab 3	60.33	10.08	29.59	100.00
Lab 4	60.37	10.05	29.51	99.93
Lab 5	60.4	10.00	29.63	100.03
Lab 6	60.32	10.12	29.62	100.06
Lab 7	60.28	10.12	29.62	100.02
Lab 8	60.33	10.08	29.58	99.99
mean	60.34	10.07	29.59	100.00
SD	0.040	0.041	0.038	0.038

Multi-element bracketing for alloys

Synthetic sample - Au 45 - Ag 845 -Pt 18 - Pd 22

Test measurements by means of a synthetic precious metal alloy using different spectral lines and measurement concentrations.

	Au			Ag			Pt			Pd		
Weight	100 mg	150 mg	300 mg	100 mg	150 mg	300 mg	100 mg	150 mg	300 mg	100 mg	150 mg	300 mg
Average [g/kg]	45,13	45,15	45,05	845,00	845,07	844,76	17,99	18,06	18,02	21,91	21,94	21,94
s [g/kg]	0,10	0,06	0,04	0,60	0,47	0,90	0,07	0,07	0,05	0,08	0,05	0,06
sr [%]	0,22 %	0,12 %	0,10 %	0,07 %	0,06 %	0,11 %	0,37 %	0,39 %	0,25 %	0,38 %	0,22 %	0,26 %
Wavelength [nm]	242	242	242	338	338	338	265	265	265	360	360	360
Measuring concentration [mg/L]	4,5	6	12	85	130	260	1,8	2,7	5,4	2,2	3,8	7,6

Multi-element bracketing for alloys

**Example for precious metal scrap
same procedure as in single element bracketing**

	Pt	Ag	Au	Pd	Ir	Rh	Ru
g/kg	207,94	32,83	81,85	73,87	47,92	11,82	18,66
mg/l	80,90	12,77	31,85	28,74	18,64	4,60	7,26
RSD	0,07	0,12	0,13	0,10	0,23	0,26	0,46

SPECTRO ICP-OES

